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Before using a product introduced onto the market by ReGen Biologies AG, the operating surgeon is asked to carefully study the following recommendations, warnings and instructions, as well as the product-specific information (technical product description, description of the surgical technique, catalogue sheet, etc.). The product-specific information contains additional data about which the operating surgeon must be aware, e.g. risks of residual material, indications or contraindications not contained in or deviating from those contained in this package insert. ReGen Biologies AG also recommends that operating surgeons attend the appropriate user-training courses. The aforementioned documents and details of the training courses may be obtained from the responsible representatives in the various countries.

### 1. Safety instructions

#### 1.1 General instructions

The CMI is designed to function as a resorbable scaffold to facilitate host meniscus tissue regeneration in patients who have an irreparable meniscus tear or loss of meniscus tissue. The CMI provides a mechanically stable, functional replacement for the deficient tissue through the implant's absorption and replacement by patient's native tissue.

The CMI may be implanted only by surgeons who are familiar with the general problems of meniscus injury and who master the product-specific surgical technique.

The CMI's are to be stored in their original packaging, unopened, and must not be damaged. Before implantation, they must be checked for any damage that may have occurred during shipping or storage. Before sterile implants are removed from the packaging, the protective packaging must be examined for possible damage as this could impair the sterility. The expiry date for the sterility of the product must be observed. When the implant is being removed from the protective packaging, the rules of asepsis must be observed.

Any complications or other effects that may occur for reasons such as incorrect indication or surgical technique, unsuitable choice of material or treatment, inappropriate use or handling of the instruments, sepsis, etc., fall under the responsibility of the operating surgeon and cannot be blamed on the manufacturer, the importers or the suppliers of products of ReGen Biologies, Inc.

The manufacturer and the distributor do not accept liability in the case of noncompliance with this package insert.

#### Important instructions:

- The CMI may only be implanted with the corresponding instruments. Instrumentation used should correspond to that described in the surgical technique.
- Instruments are always components of a system. Care should be taken to use the corresponding measuring and delivery devices of the CMI. No liability is accepted for instruments of third parties that are used by the purchaser or the user.
- The CMI's must not be altered in any way and adherence to the surgical technique must always be demonstrated.
- The general risks involved in implants are allergic reactions to the implant material used and improper positioning of the implant.
- The use of the implants for other purposes and non-approved indications is forbidden.
- Any additional instructions (e.g. instruction labels on the packaging) are to be followed.

#### Warnings:

- A CMI that has been removed from the sterile package must not be re-used.
- CMI's that are contaminated, unsterile, damaged, torn or have been improperly handled or altered without authorization must not be implanted under any circumstance.
- CMI's may tear or have impaired functionality if they are subjected to excessive loading, have been damaged, implanted incorrectly or handled improperly.
- This device is derived from animal tissue (including bovine) and although care has been exercised in its manufacture to minimize the human pathogen content, there is always a risk of the presence of residual infective agents. Eliminating this risk is a major goal of the strict manufacturing, testing and regulatory approval processes.

### 1.2 Cleaning and sterilization

The CMI is gamma irradiated.

Before the implant is used, the expiry date must be noted and the protective packaging and the sterile packaging must be checked for damage.

### 1.3 Resterilization

The CMI implant cannot be resterilized.

### 1.4 Indications

- Prior loss of meniscus tissue.
- Irreparable meniscus tears requiring partial meniscectomy.
- Either traumatic or chronic posttraumatic meniscus tear.
- Meniscus damage requiring greater than 25% removal.
- Intact anterior and posterior attachments of the meniscus.
- Intact rim over the entire circumference of the involved medial meniscus.
- ACL deficiencies corrected within 12 weeks of CMI surgery.
- Patients willing to follow post-operative rehabilitation program. Patients should be capable of understanding and following the doctor's instructions.

### 1.5 Contraindications

- Concomitant PCL insufficiency of the involved knee.
- Diagnosis of grade IV degenerative cartilage disease in the affected joint.
- Uncorrected malformations, axial malalignment in the lower extremity.
- Documented allergy to collagen of animal origin.
- Documented allergy to chondroitin sulfate of animal origin.
- Systemic or local infection.
- History of anaphylactoid reaction.
- Systemic administration of any type of corticosteroid, antineoplastics, immunosuppressing or immunosuppressive agents within 30 days of surgery.
- Evidence of osteonecrosis in the involved knee.
- Medical history that is positive for, but not restricted to, the following diseases:
  - Rheumatoid arthritis.
  - Relapsing polyarthritis.
  - Severe degenerative osteoarthritis.
  - Inflammatory arthritis.
- General neurological abnormalities or neurological conditions that tend to prevent the patient's mental ability or willingness to fulfill the requirements of the rehabilitation program.

### 1.6 Precautions

- The safety and effectiveness of the CMI in pregnant women and nursing mothers have not been established.
- The safety and effectiveness of the CMI in patients who are skeletally immature (less than 18 years of age or no radiographic evidence of epiphyseal closure) have not been established.
- The safety and effectiveness using all-inside meniscus fixation devices to fix the CMI to the native meniscus rim have not been established.
- It is not recommended to release the lateral collateral ligament when implanting the CMI in the lateral compartment.
- Avoid suturing the popliteal tendon when implanting the CMI into the lateral compartment.

### 1.7 Factors that can jeopardize the success of the operation

- History of drug or substance abuse.
- Evidence of osteoporosis or local bone tumors in the involved extremity.
- Systemic diseases and metabolic disorders.
- Pregnancy.
- Allergy to the implanted material, mainly collagen.
- Overweight/obesity (defined as patients whose weight exceeds 30% of the normal weight for their height or patients weighing more than 100 kg).
- Osteoporosis or osteomalacia.
- Unbalanced diet, medication abuse, smoking, alcohol consumption, use of drugs.
- Following the operation, physical activity should be limited to the rehabilitation protocol.
- Patients who are incapable of understanding and following the doctor's instructions.

### 1.8 Possible negative effects

The effects listed below are among the most frequent negative consequences of an implantation:

- Transient or permanent nerve damage due to pressure or haematoma.
- Complications associated with arthroscopic technique.
- Wound haematoma and slow wound healing.
- Pain.
- Restricted freedom of movement.
- Any complication associated with arthroscopy including, but not limited to, damage to articular cartilage, healthy meniscus tissue, neurovascular structures of the knee.
- Any complications associated with general surgery including blood loss, formation of blood clots (deep vein thrombosis), infection, neurological and/or vascular injury, pulmonary embolism and cardiac arrest, etc.
- Any negative effects associated with drugs or ancillary devices used during implantation. Refer to information provided by the manufacturer of the specific drug or device.

### 1.9 Information for the patient

The doctor must explain the risks of implantation of a CMI to the patient, including the possible effects of the factors mentioned under Point 1.7 on the success of the operation, and the possible negative effects mentioned under Point 1.8. The patient must also be informed as to what measures he/she can take in order to reduce the possible effects of these factors.

### 2. Surgical technique

We strongly recommend the following:

- The meniscus rim must be appropriately prepared: any unstable and degenerative meniscus tissue should be removed and a bleeding bed carefully prepared.
- The CMI should be sutured in place: vertical slits placed every 4.5 mm are recommended. Horizontal slits should be used in the posterior and anterior junctions. We recommend that non-resorbable sutures should be used except in the immediate area of the popliteal tendon. Avoid suturing the popliteal tendon when implanting the CMI into the lateral compartment. ReGen recommends use of Smith & Nephew Fast-Fix™ all-inside suture system limited only to fixation where suture placement is difficult.
- The surgeon must take precautions to avoid damaging the surrounding neurovascular structures.

Refer to the CMI Surgical Technique for additional details.

### 3. Implant materials

Bovine Type I collagen, chondroitin sulfate, and sodium hyaluronate.

### 4. Storage and handling of implants and instruments

- Careful handling is required to avoid damage to the implant.
- The CMI must be stored in the original packaging, unopened.
- The package containing the CMI must be stored at temperatures between 2°C and 25°C.

### 5. Pictograms

- "Follow the instructions for Use"
- "Not to be re-used"
- "To be used by... (Year, Month)"
- "Sterile" and "Sterilization by irradiation"
- "Temperature limitation"

### 6. Trademarks

©2006 ReGen® is a registered trademark and CMI™ is a trademark of ReGen Biologies, Inc. Redwood City, CA 94063, USA and Franklin Lakes, NJ 07417, USA.

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